

Easy setup, easy operation.

The Game Ready System is easy to set up, with an intuitive interface and fingertip controls for treatment customization. An optional battery pack enables treatment on the go.



FILL ICE AND WATER



SELECT PROGRAM (OPTIONAL)



ADJUST SETTINGS



START TREATMENT

Want to improve outcomes, satisfaction, and your bottom line? Please contact your Game Ready® rep or visit gameready.com.

AVANOS | GAME READY®

SUPPORT POSITIVE RECOVERY AFTER WORKPLACE INJURIES



GAME READY®

Get next-level cold and compression therapy.

The Game Ready® System combines a powerful control unit with a complete range of anatomic, dual-action wraps to bring proven cold and intermittent compression therapies together with all-new adjustability, precision, and convenience.

*"Reinforces treatment compliance and reduces the use of narcotics - it helps our employees come back faster."**

Laura Schwartz, RN, MSN, APRN, FNP, COHNS



*Based on individual experience, not a manufacturer's claim, results may vary. **Thermal images and skin temperatures taken immediately after 30 minutes of treatment, all devices applied per manufacturer's instructions under identical conditions.

1. Su EP, Perna M, Boettner F, et al. A prospective, multi-center randomized trial to evaluate the efficacy of a cryopneumatic device on total knee arthroplasty recovery. J Bone & Joint Surg. 2012;SuppA(94-B):153-6. 2. Waterman W, Walker JJ, Swains C, et al. The efficacy of combined cryotherapy and compression compared with cryotherapy alone following anterior cruciate ligament reconstruction. J Knee Surg. 2012;25(02):155-160. 3. Murgier J and Cassard X. Cryotherapy with dynamic intermittent compression for analgesia after anterior cruciate ligament reconstruction. Preliminary study. Orthop Traumatol Surg Res. 2014;100(3):309-312. 4. As with all RICE applications. 5. Leegwater NC, Jore HW, Brohet R, Nolte P. Cryocompression therapy after elective arthroplasty of the hip. Hip Int. 2012;22(05):527-533.

There are inherent risks in all medical devices. Please refer to the product labeling for Indications, Cautions, Warnings and Contraindications. Refer to www.gameready.com for product safety technical bulletins.

GAME READY®
ACTIVATE NEXT-LEVEL INJURY & POST-OP RECOVERY



Active Cryotherapy. Circulating ice water continuously and consistently cools the injury site to safely and comfortably reduce pain without narcotics.



Intermittent Compression. Gently squeezes and releases to remove swelling, stimulate the flow of blood, and conform the wrap to the body.



Dual-Action Wraps. Anatomically designed, circumferential wraps surround the injury site, covering more surface area for deeper cooling.



Simple Controls. The system makes it easy to use pre-set programs and quickly customize temperature, compression, and time settings.

A SYSTEM THAT HELPS REDUCE PAIN, SWELLING, AND THE USE OF OPIOIDS^{1,2,3,4}

Following workplace injuries or accidents, positive recoveries are an important step towards getting back to the things that matter. Earlier returns to work are better for people and better for your organization. The Game Ready® System uniquely integrates cold and compression therapies to support outcomes, all while helping to reduce the use of opioid medications.

Faster, deeper, longer-lasting cooling for proven therapeutic benefits:

- ✓ **Reduces opioid consumption** and may allow for **earlier discontinuation of opioid medications** in certain patients.^{1,2,3}
- ✓ **Decreases pain and swelling.**^{1,4}
- ✓ **Can increase post-operative range of motion.**³
- ✓ **May increase patient satisfaction** with the recovery process.^{1,5}

*“With Game Ready, we’ve seen decreases in post-operative pain and hospital stays, and quicker returns to work.”**

Sanford Kunkel, MD | Orthopedic Surgeon


Home recovery with Game Ready is made easy with a full-service, hassle-free rental program.



Game Ready Rental Program. A full-service rental program takes the hassle out of paperwork, insurance billing, patient education, and shipping logistics - so recovery at home is just as easy as at an on-site facility.

Trusted by the world’s top physicians and physical therapists, professional sports teams, and special military forces, the Game Ready System takes traditional RICE methods (Rest-Ice-Compression-Elevation) to the next level.

THERMAL IMAGERY SHOWS THAT GAME READY OUT-COOLS THE COMPETITION**



Innovative wrap design for better comfort and results.

Anatomically engineered for all major body parts, Game Ready’s patented wraps ensure better fit, greater comfort, and circumferential coverage for more effective therapy.

